
PenEvac® Warnings, Cautions & Instructions

I.C. Medical, Inc.
 15002 North 25th Drive
 Phoenix, Arizona 85023 USA
 Phone: (623) 780-0700 (USA)
 Fax: (623)780-0887 (USA)
 www.icmedical.com

Made in USA

en

Indications for Use

The PenEvac® is intended to be used as the active **Monopolar** electrode in an electrosurgery generator system and to facilitate the removal of smoke that is generated during procedures.

Maximum electrical capacity:

The PenEvac is intended to be used with an Electrical Surgical Unit. The maximum electrical capacity of this product is 4.5 kVp. **DO NOT** exceed this value under any circumstances. The PenEvac is intended to be used with an I.C. Medical smoke evacuator.

Safety Instructions

WARNINGS!

- Keep active accessories away from patient when not in use and away from flammable objects, gases, and vapors at all times
- This is a monopolar device. BEFORE USE, a dispersive electrode must be properly attached to the patient.
- Keep active accessories in a clean, dry, non-conductive safety container (holster) provided when not in use.
- After using the Cut or Coagulation function on an electrosurgical pencil (PenEvac®), the tip/electrode is hot. **DO NOT** lay any active device on the patient, especially an electrosurgical pencil (PenEvac®) that has just been used on the patient. This is to prevent the patient, physician, and/or staff from accidental burns by the hot tip/electrode.
- Electrosurgical electrodes that are activated or hot from use can cause a fire. Do not place them near or in contact with flammable materials and substances (e. g. drapes, linens, flammable gases, endotracheal tubes, etc.)
- Keep the voltage /power as low as possible to achieve the desired end effect in order to minimize the potential for capacitive coupling and inadvertent burning at high voltages.
- Check the product for damage before use, including all the insulation. If damaged, **DO NOT** use this product.
- **DO NOT USE** in MR Environment!
- **DO NOT** use this device for the suction of liquid
- **DO NOT USE** in the presence of flammable anesthetics, oxidizing gases (such as nitrous oxide (N₂O)

and oxygen) or any potential source of ignition!

- Protect this product from any form of mechanical damage! **DO NOT** throw! **DO NOT** use force! If there is a cable, and unless stated otherwise, **DO NOT** kink it or wrap it around the product!
- **DO NOT** kink the tube or wrap it around the product.
- **DO NOT USE** the telescope for lifting tissue, cracks may appear and electrical integrity may be compromised.

CAUTIONS:

- Federal (U.S.A.) Law restricts this device to sale by or on the order of a licensed physician.
- For all PenEvac® Accessories: Read instructions for PenEvac®, Electrosurgery Unit, and Smoke Evacuator before Use.
- **DO NOT** operate the electrosurgical pencil if air does not flow through the smoke evacuator tube. Refer to smoke evacuator instructions
- **DO NOT USE** a scratch pad or other abrasive cleaner to remove eschar. This may damage the PTFE coating.
- If the electrode or coating is damaged, discard the electrode.
- Activate the electrodes when in contact, or in close proximity to target tissue to avoid the possibility of unintended tissue damage.
- Electrosurgical cables should be positioned to minimize contact with the patient and avoid contact with other leads to avoid adversely influencing the operation of other electronic equipment.
- **DO NOT USE** in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional

(e.g. cardiologist), to avoid an interference with the action of the electronic implant. The implant may be damaged, which would imply a possible risk for the patient.

- **DO NOT** use for fluid removal, this device was not designed for such application. Aspirate fluid from the area before activating the instrument. Direct contact or near an active electrode with conductive fluids (e.g., blood or saline) may carry electrical current or heat away from target tissue, which may cause burns to the patient.
- Localized burns to the patient or physician may result from electric currents carried through conductive objects. Current may be generated in conductive objects by direct contact with the active electrode, or by the active accessory being in close proximity to the conductive object.
- When the visualization may be impaired, be alert, inadvertent activation or movement of the activated electrode, may result in injury to the patient.

INSTRUCTIONS:

Sterile product for single use

Single-use device. Properly discard after use. **DO NOT** re-sterilize or reuse. Reuse may result in patient injury due to contamination or insulation failure.

The sterilization method used for this product is ETO. For details on sterilization process, contact I.C. Medical, Inc.

Non-Sterile product

The "Non-Sterile" product is indicated on the packaging by Non-Sterile symbol.

I.C. Medical intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient.

This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

Warning: Malfunctions:

If the device malfunctions, immediately discontinue use.

Expiration date:

The use-by-date is indicated on the packaging. **DO NOT** use this product if the use-by-date has been exceeded!

Packaging:

- Use aseptic techniques when removing the product from the packaging!
- Device sterilized by EO gas.
- Use this product only if the packaging has not been opened or damaged. If the packaging is open or damaged, the product is no longer sterile. In this case throw the product away! **DO NOT re-sterilize.**

Operating Instructions:

- Check the product for damage before use, including

all the insulation. If damaged **DO NOT** use this product

- Insert electrical connector to **Monopolar** connector on ESU generator.
- Connect the smoke tube to the smoke evacuator.
- On the Telescopic PenEvac® loosen the locking ring ¼ turn counterclockwise, while looking into the nozzle, adjust to desired length and tighten the locking ring on the telescopic electrode.
- Note: When the locking ring is secure, the electrode can still be rotated to the desired orientation without loosening the ring. To change the length of the telescoping tip, you will still have to loosen the ring.

Cleaning, disinfection, sterilization

DO NOT clean, **DO NOT** disinfect, **DO NOT** re-sterilize or reuse!

PenEvac® Electrodes (All Styles) Instructions:

- Use aseptic technique. Open package.
- Remove electrode from the PenEvac®
- Insert new electrode in the socket inside the PenEvac®. Make sure that the electrode is securely inserted in the socket.
- Adjust length of telescope and tighten locking ring.

Smoke Tube Instructions:

- Use aseptic technique (if product is labeled "Sterile"). Open Package.
- Connect tubing to smoke evacuator filter and to the handpiece connector.
- If using 7/8" Smoke Tubing, connect to smoke evacuator filter. Connect other end to smoke evacuator accessory, using any required adapters.


- Single use devices. **DO NOT** clean, disinfect, re-sterilize or reuse.

All Items

- Replace and dispose all items labeled "Single Use" in accordance with your institutional biohazard protocol.
- **Do Not clean, re-sterilize or reuse.** Reuse may result in patient injury due to contamination or dielectric failure.
- Transport and storage conditions: Temperature limit -18°C to 60°C (0°F to 140°F); Humidity limitation 30% to 85%.

Serious Adverse Events:

Any serious adverse event or incident that occurs in relation to the device or accessory should be reported to the manufacturer, I.C. Medical, Inc., at complaints@icmedical.com and to the FDA. In addition, European customers should also report to the Authorized Representative at the address listed on the label or IFU and to the competent authority in the member state.

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