I.C. Medical, Inc.
 Manufacturer 15002 North 25th Drive
 Phoenix, Arizona 85023 USA
 Phone: (623) 780-0700 (USA)
 Fax: (623)780-0887 (USA)

Made in USA

EC REP MDSS GmbH
 Schiffgraben 41
 30175 Hannover,
 Germany

en

**Smoke Evacuator
 Accessory Warnings,
 Cautions & Instructions**

ESU SHROUD

MD **CE** **2797** **Caution** **Do not re-use**

140°F 60°C
 0°F -18°C Temperature limit 30% Humidity limitation 85%

Non-Sterile **Rx Only**

Not manufactured using natural rubber latex

Indications for Use

The ESU Shroud, slips over standard electrocautery (ESU) hand switching pencil and is used to evacuate smoke and other airborne debris that is created when the ESU pencil is in use.

GENERAL ESU

WARNINGS:

- Keep active accessories away from patient when not in use and away from flammable objects, gases, and vapors at all times.
- Keep active accessories in the safety container (holster), when not in use.

• After using the Cut or Coagulation function on an electrocautery pencil, the tip/electrode is hot. Do not lay any active device on the patient, especially an electrocautery pencil that has just been used on the patient. This is to prevent the patient, physician, and/or staff from accidental burns by the hot tip/electrode.

• Keep the voltage/power as low as possible to achieve the desired end effect in order to minimize the potential for capacitive coupling and inadvertent burning at high voltages.

• Check the product for damage before use, including all the insulation. If damaged, DO NOT use this product.

• Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! If there is a cable, and unless stated otherwise, DO NOT kink it or wrap it around the product!

• DO NOT use this device for the suction of liquids

• Use only under the direction of a licensed physician.

• The ESU Shroud may be used with an I.C. Medical Smoke Evacuator and any standard electrocautery (ESU) hand switching pencil having an approximate length of 6in to 6.75in and outside diameter of 0.4in to 0.5in.

WARNINGS:

- Check the product for damage before use. If damaged, DO NOT use this product.
- Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force!
- Use only under the direction of a licensed physician.

CAUTIONS:

- Read instructions for Electrocautery Unit (ESU)

and Operating and Installation Manual for Smoke Evacuator.

• **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

INSTRUCTIONS:

• Place electrocautery pencil into shroud sleeve. Make sure it is all the way forward. Pencil should be snug in shroud and pencil blade tip should protrude approximately 1/3" to 3/4" beyond the end of the shroud tip.

• The shroud is supplied with a standard-length tip suitable for most uses.

• For surgeon's comfort, fasten tubing (about 3 feet from pencil end) to a drape near the surgical site. Use towel clamp, or similar device. Do not occlude tubing.

• Attach the hose to the input filter of the I.C. Medical smoke evacuator.

• Connecting the electrocautery, pencil cord to the shroud hose will make it easier to move the shroud and pencil during surgery.

• Air should flow through the shroud when pencil is active. If not, check for obstructions and see smoke evacuator operator manual.

• Replace and dispose all items labeled "Single Use" in accordance with your institutional protocol. **Do Not clean, re-sterilize or reuse.**

• Transport and storage conditions: Temperature limits -18°C to 60°C (0°F to 140°F); Humidity limitation 30% to 85%

Sterile product

2 **STERILE EO**

Do not re-sterilize Do not use if package is damaged Sterilized using ethylene oxide.

The sterilization method used is ETO. For details on

sterilization process, contact I.C. Medical, Inc.

PACKAGING:

• Use aseptic techniques when removing the product from the packaging.

• Use this product only if the packaging has not been opened or damaged. In this case throw the product away

Non-Sterile product

The product "Non-Sterile" is indicated on the packaging by

Non-Sterile symbol.

I.C. Medical intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient.

WARNINGS:

• Check the product for damage before use. If damaged, DO NOT use this product.

• Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force!

• Use only under the direction of a licensed physician.

CAUTIONS:

• Read Operating and Installation Manual for Smoke Evacuator.

• Before installing or removing any filter, be sure that the system is turned off and the electrical plug has been disconnected.

• Failure to install the SAFEGUARD BLUE® Hydrophobic ULPA Filter with Built-in Fluid Trap could result in contamination and possible electrical damage of the I. C. Medical Smoke Evacuator.

• **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

MD **CE** **Caution** **Do not re-use**

158°F 70°C
 40°F -40°C Temperature limit 100% Humidity limitation 10%

Non-Sterile **Rx Only**

Not manufactured using natural rubber latex

Indications for Use

The filter is used as an attachment for all I.C. Medical Smoke Evacuator models, to remove, filtrate the smoke plume and to collect in its built-in fluid trap, any incidental fluid, vapors, oil mist, hydrocarbon vapors, fat vapors, moisture droplets, odors, organic solvent vapors, and toxic vapors from the surgical site.

WARNINGS:

• Check the product for damage before use. If damaged, DO NOT use this product.

• Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force!

• Use only under the direction of a licensed physician.

CAUTIONS:

• Read Operating and Installation Manual for Smoke Evacuator.

• Before installing or removing any filter, be sure that the system is turned off and the electrical plug has been disconnected.

• Failure to install the SAFEGUARD BLUE® Hydrophobic ULPA Filter with Built-in Fluid Trap could result in contamination and possible electrical damage of the I. C. Medical Smoke Evacuator.

• **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

INSTRUCTIONS:

• The SAFEGUARD BLUE® Hydrophobic ULPA Filter with Built-in Fluid Trap should only be used with I.C. Medical Smoke Evacuators (all models).

• Attach the filter assembly to the smoke evacuator and attach the accessory tubing to the filter assembly.

• Replace the filter assembly when "Change Filter" indicator comes on.

• Replace the filter assembly when desired levels of gas flow cannot be achieved.

• Replace the filter assembly after any procedures when moisture or fluids are present in the water/vapor trap.

• Replace cap on end of filter when not in use.

• Dispose of filter according to your institution's policy for biologically contaminated waste disposal.

• This is a NON-STERILE product. It may be used until there is an indication to replace the filter (listed above). After that time do not attempt to clean, sterilize, or re-use this filter.

• Transport and storage conditions: Temperature limit -40°C to 70°C (-40°F to 158°F); Humidity limitations 10% to 100%.

PACKAGING:

• Use this product only if the packaging has not been opened or damaged. In this case throw the product away!

**INTRA-ABDOMINAL
 PLUME ELIMINATOR
 TUBING SETS**

MD **CE** **Caution** **Do not re-use**

140°F 60°C
 0°F -18°C Temperature limit 85% Humidity limitation 30%

Non-Sterile **Rx Only**

Not manufactured using natural rubber latex

Indications for Use

It is used to eliminate smoke and any airborne debris that are generated during laparoscopic procedures.

WARNINGS:

• Check the product for damage before use. If damaged, DO NOT use this product.

• Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! DO NOT kink it or wrap it around the product!

• Use only under the direction of a licensed physician

CAUTIONS:

• Read Operating and Installation Manual for Smoke Evacuator.

• **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

INSTRUCTIONS:

• The product connects to any model of an I.C. Medical Smoke Evacuator at one end and to laparoscopic trocar at the other end.

• Firmly attach adapters to desired connectors.

• Replace and dispose all items labeled "Single Use" in accordance with your institutional protocol. **Do Not clean, re-sterilize or reuse.**

• Transport and storage conditions: Temperature limit -18°C to 60°C (0°F to 140°F); Humidity limitation 30% to 85%.

Sterile product

CE **2** **STERILE EO**

2797 Do not re-sterilize Do not use if package is damaged Sterilized using ethylene oxide.

The sterilization method used is ETO. For details on sterilization process, contact I.C. Medical, Inc.

Non-Sterile product

CE **Non-Sterile**

The product "Non-Sterile" is indicated on the packaging by Non-Sterile symbol.

I.C. Medical intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient.

This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death

PACKAGING:

• Use aseptic techniques when removing the product from the packaging.

• Use this product only if the packaging has not been opened or damaged. In this case throw the product away!

SPECULUM TUBING SETS

MD **CE** **Caution** **Do not re-use**

140°F 60°C
 0°F -18°C Temperature limit 85% Humidity limitation 30%

Non-Sterile **Rx Only**

Not manufactured using natural rubber latex

Indications for Use

Smoke Accessories are intended to evacuate smoke plume produced during surgical procedures.

WARNINGS:

• Check the product for damage before use. If damaged, DO NOT use this product.

• Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! DO NOT kink it or wrap it around the product!

CAUTIONS:

• Read Smoke Evacuator Accessory Warnings, Cautions & Instructions; Electrocautery Unit (ESU); and Operating and Installation Manual for Smoke Evacuator.

• **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

INSTRUCTIONS:

• Attach one end of speculum tubing to the speculum port, and opposite end of the tubing to the SAFEGUARD BLUE® Hydrophobic ULPA Filter with Built-in Fluid Trap on the I.C. Medical smoke evacuator.


• Replace and dispose all items labeled "Single Use" in accordance with your

Smoke Evacuator Accessory Warnings, Cautions & Instructions

institutional protocol. **Do Not clean, re-sterilize or reuse.**

- Transport and storage conditions: Temperature limit -18°C to 60°C (0°F to 140°F); Humidity limitation 30%-85%.

Sterile product

The sterilization method used is ETO. For details on sterilization process, contact I.C. Medical, Inc.

Non-Sterile product

The product "Non-Sterile" is indicated on the packaging by Non-Sterile symbol.


I.C. Medical intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient. This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

PACKAGING:

- Use aseptic techniques when removing the product from the packaging.
- Use this product only if the packaging has not been

opened or damaged. In this case throw the product away!

SMOKE EVACUATOR WANDS

Indications for Use

Smoke Accessories are intended to evacuate smoke plume produced during surgical procedures.

WARNINGS:

- Check the product for damage before use. If damaged, DO NOT use this product.
- Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force. DO NOT kink it or wrap it around the product!
- Use only under the direction of a licensed physician.

CAUTIONS:

- Read Operating and Installation Manual for Smoke Evacuator.
- **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.


INSTRUCTIONS:

- Attach wand hose to the input filter of the smoke evacuator.
- Place wand close to source of smoke.
- Replace and dispose all items labeled "Single Use" in

accordance with your institutional protocol. **Do Not clean, re-sterilize or reuse.**

- Transport and storage conditions: Temperature limit -18°C to 60°C (0°F to 140°F); Humidity limitation 30% to 85%.

Sterile product

The sterilization method used is ETO. For details on sterilization process, contact I.C. Medical, Inc.

Non-Sterile product

The product "Non-Sterile" is indicated on the packaging by Non-Sterile symbol.

I.C. Medical intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient. This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.


I.C. Medical intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient. This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

PACKAGING:

- Use aseptic techniques when removing the product from the packaging.

Use this product only if the packaging has not been opened or damaged. In this case throw the product away!

ULPA IN-LINE SUCTION FILTERS WITH or WITHOUT TUBING SETS

Indications for Use

The filter is used as an attachment to smoke accessories or smoke evacuator devices to filter surgical smoke, odors, organic solvent vapors, toxic vapors, noxious odors produced during surgical procedures.

WARNINGS:

- Check the product for damage before use. If damaged, DO NOT use this product.
- Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! DO NOT kink it or wrap it around the product!
- Use only under the direction of a licensed physician.

CAUTIONS:

- Read Operating and Installation Manual for Smoke Evacuator.
- **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

WARNINGS:

- Check the product for damage before use. If damaged, DO NOT use this product.


INSTRUCTIONS:

- Attach to smoke accessories, suction walls, smoke tubing or smoke evacuator devices.
- Replace the filter assembly when discoloration shows on filter.
- Dispose of filter according to your institution's policy for biologically contaminated waste disposal.
- Transport and storage conditions: Temperature limit -18°C to 60°C (0°F to 140°F); Humidity limitation 30% to 85%.
- This is a NON-STERILE product. It may be used until there is an indication to replace the filter (listed above). After that time do not attempt to clean, sterilize, or re-use this filter.

PACKAGING:

- Use this product only if the packaging has not been opened or damaged. In this case throw the product away!

SMOKE TUBES

Indications for Use

Smoke Accessories are intended to evacuate smoke plume produced during surgical procedures.

WARNINGS:

- Check the product for damage before use. If damaged, DO NOT use this product.

- Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! DO NOT kink it or wrap it around the product!
- Use only under the direction of a licensed physician.


CAUTIONS:

- Read Operating and Installation Manual for Smoke Evacuator.
- **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

INSTRUCTIONS:

- Connect tubing to smoke evacuator's filter and to the handpiece connector.
- If using 7/8" Smoke Tubing, connect to smoke evacuator filter. Connect other end to smoke evacuator accessory, using any required adapters.
- Replace and dispose all items labeled "Single Use" in accordance with your institutional protocol. **Do Not clean, re-sterilize or reuse.**
- Transport and storage conditions: Temperature limit -18°C to 60°C (0°F to 140°F); Humidity limitation 30% to 85%.

Sterile product

The sterilization method used is ETO. For details on sterilization process, contact I.C. Medical, Inc.

Non-Sterile product

The product "Non-Sterile" is indicated on the packaging by Non-Sterile symbol. I.C. Medical intends that non-sterile products for use in

sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient. This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

PACKAGING:

- Use aseptic techniques when removing the product from the packaging.
- Use this product only if the packaging has not been opened or damaged. In this case throw the product away!

SERIOUS ADVERSE EVENTS:

Any serious adverse event or incident that occurs in relation to the device or accessory should be reported to the manufacturer, I.C. Medical, Inc., at complaints@icmedical.com and to the FDA. In addition, European customers should also report to the Authorized Representative at the address listed on the label or IFU and to the competent authority in the member state.

I.C. Medical intends that non-sterile products for use in

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