
I. C. MEDICAL, INC.

LASER SENSOR INSTALLATION INSTRUCTIONS

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MDSS-UK RP Ltd.
6 Wilmslow Road,
Rusholme,
Manchester M14 5TP,
United Kingdom

Made in the USA

INDICATIONS FOR USE:

The Laser Sensor is used as an attachment to a laser foot switch to turn the smoke evacuator ON and OFF.

TECHNICAL SPECIFICATION:

I.C. Medical, Inc. recommends the use of Laser Sensor with the Crystal Vision® Smoke Evacuators-all models. Laser Sensor can be installed on the housing of a laser foot switch and connected to the Crystal Vision Smoke Evacuator to enable activation during surgical procedures.

TRANSPORT AND STORAGE CONDITIONS:

Temperature limits: -40°C to 70°C (-40°F to 158°F)

Relative Humidity limitation: 10% to 100%.

Warranty period is 2 months from date of sale.

Without limitation, the warranty does not cover damage caused by customer misuse or abuse of a product.

LASER SENSOR INSTALLATION

If you are going to use your CRYSTAL VISION® with a laser, proceed with the following:

1. Attach the Sensor Cable (Figure 1) to the LASER connector on the CRYSTAL VISION® and to the Sensor Connector Box.
2. The Sensor Assembly has three parts: the Sensor Connector Box, the Transmit Sensor, and the Receive Sensor. Both Sensors have double-backed tape on one side and an infrared lens on the side opposite the tape. The Receive Sensor has a red indicator lamp that will light when the sensor assembly is plugged into the operating CRYSTAL VISION®. (NOTE: When the red light goes off, CRYSTAL VISION® will start to operate).
Test the Sensor Assembly for proper operation:
 - Plug the Sensor Cable into the Sensor Connector Box and the LASER connector of the CRYSTAL VISION®.
 - Plug the CRYSTAL VISION® into an AC outlet and turn on the Power Switch on the Control Box back panel and front panel.
 - Align the Transmit Sensor Lens and the Receive Sensor Lens until the red light goes off and the CRYSTAL VISION® starts.
 - Move the Sensor until the red light comes on. The pump in the CRYSTAL VISION® will stop operating sometime within 30 seconds after the light goes out. (NOTE: the actual amount of time that is required for the pump to cease operation is determined by the TIME adjustment on the CRYSTAL VISION® front panel).
3. Place the foot switch for the laser on an easily accessible work surface (Figures 2 & 3).
4. Inside the foot switch housing (Figure 2), position the Transmit Sensor and the Receive Sensor on opposite sides of the foot switch. Do not position the sensors on the sides of the foot switch. **WITHOUT REMOVING THE PROTECTIVE COVERING ON THE TAPE**, position them so that the red light comes on.
5. Keep the Sensors in the same position and press the Laser Foot Pedal down. The red light should go out. If it does not, reposition the Sensors until the red light goes out when the Laser Foot Pedal is depressed and it remains on when the pedal is not depressed.
THE LIGHT SHOULD GO OUT FOR THE SLIGHTEST MOVEMENT OF THE FOOT PEDAL. If it does not, move both sensors higher up the wall of the protective housing of the laser pedal. Carefully mark the location of both sensors.
6. Remove the protective backing from one Sensor and place it in the correct position on the side of the Laser Foot Switch Assembly. (NOTE: It is usually very helpful to only LIGHTLY position the sensors at first and only after you are certain that they are in the exact position press them firmly into place.)
WARNING- Repositioning either sensor after it has been firmly set in place can easily destroy it
7. Repeat the process for the other Sensor.

8. Move the Sensor Connector Box to a convenient location inside the Laser Foot Switch Assembly. Be sure that the cables from the Sensors to the Sensor Connector Box do not interfere with the operation of the foot switch or with the surgeon's foot. Carefully mark this location
9. Remove the protective covering from the tape on the Sensor Connector Box and attach it to the previously marked location inside the Laser Foot Switch Assembly.
10. Keep active accessories away from patient when not in use and away from flammable objects, gases, and vapors always.
11. Keep active accessories in the safety container (holster), when not in use.
12. After using the Cut or Coagulation function on an electrosurgical pencil, the tip/electrode is hot. Do not lay any active device on the patient, especially an electrosurgical pencil that has just been used on the patient. This is to prevent the patient, physician, and/or staff from accidental burns by the hot tip/electrode.
13. Keep the voltage/power as low as possible to achieve the desired end effect in order to minimize the potential for capacitive coupling and inadvertent burning at high voltages.
14. Check the product for damage before use, including all the insulation. If damaged, DO NOT use this product.
15. Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! If there is a cable, and unless stated otherwise, DO NOT kink it or wrap it around the product!
16. Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! If there is a cable, and unless stated otherwise, DO NOT kink it or wrap it around the product!
17. "Studies have shown that smoke generated during electrosurgical procedures may be harmful to surgical personnel. These studies recommend using a surgical mask and adequate ventilation of the smoke using a smoke evacuator or other means" (ANSI/AAMI HF-18: 1993).
18. Use only under the direction of a licensed physician.
19. The Smoke Evac and Accessory products are intended for use only with other IC Medical products. Using any other Smoke Evac or accessory not supplied by IC Medical

SERIOUS ADVERSE EVENTS:

Any serious adverse event or incident that occurs in relation to the device or accessory should be reported to the manufacturer, I.C. Medical, Inc., at complaints@icmedical.com and to the FDA. In addition, European customers should also report to the Authorized Representative at the address listed on the label or IFU and to the competent authority in the member state.

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