IC. MEDICAL, INC.

ESU SHROUD, WARNINGS, CAUTIONS & INSTRUCTIONS

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Australian Sponsor:
I.C. Medical Australia Pty Ltd
Level 1, Suite 3, 130 Bundall Road
Bundall, QLD 4217,
Australia

Made in the USA

Label: ICM-001-0219 Rev. S Date: 2025-10

Indications for Use

The ESU Shroud, slips over standard electrosurgery (ESU) hand switching pencil and is used to evacuate smoke and other airborne debris that is created when the ESU pencil is in use.

Technical Specification:

I.C. Medical, Inc. recommends the use of the ESU Shroud with the Crystal Vision® Smoke Evacuators-all models, and other I.C. Medical products.

The ESU Shroud product family, when used as a capture device positioned at a 45-degree angle and approximately 1 cm from the smoke source, with ø10 mm transfer tubing, and without interfering with the visibility of the surgical site, meets the ISO 16571:2024 standard for plume removal efficiency when paired with the Crystal Vision® Smoke Evacuator equipped with the SAFEGUARD BLUE® Hydrophobic ULPA Filter with a built-in fluid trap, operating at flow settings of ≥75 LPM.

Transport and Storage Conditions:

Temperature limits: -18°C to 60°C (0°F to 140°F);

Relative Humidity limitation: 30% to 85%.

GENERAL ESU WARNINGS:

1. Keep active accessories away from patient when not in use and away from flammable objects, gases, and vapors at all times.
2. Keep active accessories in the safety container (holster), when not in use.
3. After using the Cut or Coagulation function on an electrosurgical pencil, the tip/electrode is hot. Do not lay any active device on the patient, especially an electrosurgical pencil that has just been used on the patient. This is to prevent the patient, physician, and/or staff from accidental burns by the hot tip/electrode.
4. Keep the voltage/power as low as possible to achieve the desired end effect in order to minimize the potential for capacitive coupling and inadvertent burning at high voltages.
5. Check the product for damage before use, including all the insulation. If damaged, DO NOT use this product.
6. Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! If there is a cable, and unless stated otherwise, DO NOT kink it or wrap it around the product!
7. DO NOT use this device for the suction of liquids
8. Use only under the direction of a licensed physician.
9. The ESU Shroud may be used with an I.C. Medical Smoke Evacuator and any standard electrosurgery (ESU) hand switching pencil having an approximate length of 6in to 6.75in and outside diameter of 0.4in to 0.5in.

WARNINGS:

1. Check the product for damage before use. If damaged, DO NOT use this product.
2. Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force!
3. Use only under the direction of a licensed physician.

CAUTIONS:

1. Read instructions for Electrosurgery Unit (ESU) and Operating and Installation Manual for Smoke Evacuator.
2. **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

INSTRUCTIONS:

1. Place electrocautery pencil into shroud sleeve. Make sure it is all the way forward. Pencil should be snug in shroud and pencil blade tip should protrude approximately 1/3" to ¾" beyond the end of the shroud tip.
2. The shroud is supplied with a standard-length tip suitable for most uses.
3. For surgeon's comfort, fasten tubing (about 3 feet from pencil end) to a drape near the surgical site. Use towel clamp, or similar device. Do not occlude tubing.
4. Attach the hose to the input filter of the I.C. Medical smoke evacuator.
5. Connecting the electrocautery, pencil cord to the shroud hose will make it easier to move the shroud and pencil during surgery.
6. Air should flow through the shroud when pencil is active. If not, check for obstructions and see smoke evacuator operator manual.
7. Replace and dispose all items labeled "Single Use" in accordance with your institutional protocol. **Do Not clean, re-sterilize or reuse.**

Sterile product:

The sterilization method used is ETO. For details on sterilization process, contact I.C. Medical, Inc.

PACKAGING:

1. Use aseptic techniques when removing the product from the packaging.
2. Use this product only if the packaging has not been opened or damaged. In this case throw the product away

Non-Sterile product:

The product "Non-Sterile" is indicated on the packaging by Non-Sterile symbol.

I.C. Medical intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient. This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

SERIOUS ADVERSE EVENTS:

Any serious adverse event or incident that occurs in relation to the device or accessory should be reported to the manufacturer, I.C. Medical, Inc., at complaints@icmedical.com and to the FDA. In addition, European customers should also report to the Authorized Representative at the address listed on the label or IFU and to the competent authority in the member state.

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