# I. C. MEDICAL, INC.

## ULPA IN-LINE SUCTION FILTER, WARNINGS, CAUTIONS & INSTRUCTIONS

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**Australian Sponsor:**  
**I.C. Medical Australia Pty Ltd**  
Level 1, Suite 3, 130 Bundall Road  
Bundall, QLD 4217,  
Australia

## **Made in the USA**

en

## ULPA IN-LINE SUCTION FILTERS WITH or WITHOUT TUBING SETS

### Indications for Use:

The filter is used as an attachment to smoke accessories or smoke evacuator devices to filter surgical smoke, odors, organic solvent vapors, toxic vapors, noxious odors produced during surgical procedures.

### Technical Specification:

I.C. Medical, Inc. recommends the use of the ULPA Inline Suction Filters With or Without Tubing Sets with the Crystal Vision® Smoke Evacuators-all models, and other I.C. Medical products.

### Transport and Storage Conditions:

Temperature limits: -18°C to 60°C (0°F to 140°F);

Relative Humidity limitation: 30% to 85%.

### WARNINGS:

1. Check the product for damage before use. If damaged, DO NOT use this product.
2. Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! DO NOT kink it or wrap it around the product!
3. Use only under the direction of a licensed physician.

### CAUTIONS:

1. Read Operating and Installation Manual for Smoke Evacuator.
2. **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

### INSTRUCTIONS:

1. Attach to smoke accessories, suction walls, smoke tubing or smoke evacuator devices.
2. Replace the filter assembly when discoloration shows on filter.
3. Dispose of filter according to your institution's policy for biologically contaminated waste disposal.
4. This is a NON-STERILE product. It may be used until there is an indication to replace the filter (listed above). After that time do not attempt to clean, sterilize, or re-use this filter.

### PACKAGING:

Use this product only if the packaging has not been opened or damaged. In this case throw the product away!

### SERIOUS ADVERSE EVENTS:

Any serious adverse event or incident that occurs in relation to the device or accessory should be reported to the manufacturer, I.C. Medical, Inc., at [complaints@icmedical.com](mailto:complaints@icmedical.com) and to the FDA. In addition, European customers should also report to the Authorized Representative at the address listed on the label or IFU and to the competent authority in the member state.

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