
I. C. MEDICAL, INC.

SPECULUM TUBING SET, WARNINGS, CAUTIONS & INSTRUCTIONS

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Australian Sponsor:
I.C. Medical Australia Pty Ltd
Level 1, Suite 3, 130 Bundall Road
Bundall, QLD 4217,
Australia

Made in the USA

Label: ICM-001-1739 Rev. A Date: 2025-10

SPECULUM TUBING SETS

Indications for Use:

Smoke Accessories are intended to evacuate smoke plume produced during surgical procedures.

Technical Specification:

I.C. Medical, Inc. recommends the use of the Speculum Tubing Sets with the Crystal Vision® Smoke Evacuators—all models, and other I.C. Medical products. The Speculum Tubing Set, when used as a capture device connects directly to the speculum with ø10 mm transfer tubing, and without interfering with the visibility of the surgical site, meets the ISO 16571:2024 standard for plume removal efficiency when paired with the Crystal Vision® Smoke Evacuator equipped with the SAFEGUARD BLUE® Hydrophobic ULPA Filter with a built-in fluid trap, operating at flow settings of ≥75 LPM in OPEN mode.

Transport and Storage Conditions:

Temperature limits: -18°C to 60°C (0°F to 140°F);

Relative Humidity limitation: 30% to 85%.

WARNINGS:

1. Check the product for damage before use. If damaged, DO NOT use this product.
2. Protect this product from any form of mechanical damage! DO NOT throw! DO NOT use force! DO NOT kink it or wrap it around the product!

CAUTIONS:

1. Read Smoke Evacuator Accessory Warnings, Cautions & Instructions; Electrosurgery Unit (ESU); and Operating and Installation Manual for Smoke Evacuator.
2. **RX ONLY** Federal Law restricts this device to sale by or on the order of a licensed physician.

INSTRUCTIONS:

1. Attach one end of speculum tubing to the speculum port, and opposite end of the tubing to the SAFEGUARD BLUE® Hydrophobic ULPA Filter with Built-in Fluid Trap on the I.C. Medical smoke evacuator.
2. Replace and dispose all items labeled “Single Use” in accordance with your institutional protocol. **Do Not clean, re-sterilize or reuse.**

Sterile product:

The sterilization method used is ETO. For details on sterilization process, contact I.C. Medical, Inc.

Non-Sterile product:

The product “Non-Sterile” is indicated on the packaging by Non-Sterile symbol.

I.C. Medical intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by I.C. Medical and may damage the product which could cause a device malfunction and/or injury to the patient. This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

PACKAGING:

1. Use aseptic techniques when removing the product from the packaging.
2. Use this product only if the packaging has not been opened or damaged. In this case throw the product away!

SERIOUS ADVERSE EVENTS:

Any serious adverse event or incident that occurs in relation to the device or accessory should be reported to the manufacturer, I.C. Medical, Inc., at complaints@icmedical.com and to the FDA. In addition, European customers should also report to the Authorized Representative at the address listed on the label or IFU and to the competent authority in the member state.

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